

Mannose Receptor/CD206 Rabbit mAb


ABZ3044

Product Overview

Name	Mannose Receptor/CD206 Rabbit mAb
Accession(Primary)	P22897
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 166 kDa; Observed MW: 166 kDa
Recommended Dilution	WB-1:5000-1:50000

Target information (P22897)

Gene Name	MRC1
About Gene Symbol	MRC1 (Mannose Receptor C-Type 1) is a membrane receptor involved in adaptive immunity. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Mannose Receptor in Jurkat, C6, HeLa lysates using Mannose Receptor antibody.