

CSF1R/CD115 Receptor Rabbit mAb

ABZ3046

Product Overview

Name	CSF1R/CD115 Receptor Rabbit mAb
Accession(Primary)	P09581
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse, Rat
Calculated MW	Calculated MW: 109 kDa; Observed MW: 109 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:100

Target information (P09581)

Gene Name Csf1R

About Gene Symbol



