


Mesothelin Rabbit mAb

ABZ3047

Product Overview

Name	Mesothelin Rabbit mAb
Accession(Primary)	Q61468
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse
Calculated MW	Calculated MW: 70 kDa; Observed MW: 40,70 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20

Target information (Q61468)

Gene Name Msln**About Gene Symbol**

Western blot analysis of Mesothelin in 3T3, mouse lung lysates using Mesothelin antibody.

