

c-Met Rabbit mAb


ABZ3048

Product Overview


Name	c-Met Rabbit mAb
Accession(Primary)	P08581
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 156 kDa; Observed MW: 170 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:50-1:100

Target information (P08581)


Gene Name	MET
About Gene Symbol	



Western blot analysis of Met (c-Met) in K562, C6, 3T3 lysates using c-Met antibody.



Immunocytochemistry analysis of Met (c-Met)(green) in Hela using Met (c-Met) antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Met (c-Met) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.