

MGST3 Rabbit mAb

ABZ3050


Product Overview

Name	MGST3 Rabbit mAb
Accession(Primary)	O14880
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 17 kDa; Observed MW: 17 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; IP-1:20


Target information (O14880)

Gene Name MGST3

About Gene Symbol



Western blot analysis of MGST3 in K562, HeLa lysates using MGST3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using MGST3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.