

MIB1 Rabbit mAb


ABZ3051

Product Overview

Name	MIB1 Rabbit mAb
Accession(Primary)	Q86YT6
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 110 kDa; Observed MW: 110 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:500; ICC/IF-1:100-1:200

Target information (Q86YT6)

Gene Name	MIB1
About Gene Symbol	MIB1 (MIB E3 Ubiquitin Protein Ligase 1) is a ubiquitin E3 ligase involved in ubiquitin-proteasome pathway. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Mib1/Mindbomb in K562, Hela lysates using MIB1 antibody.