

MIF Rabbit mAb


ABZ3052

Product Overview


Name	MIF Rabbit mAb
Accession(Primary)	P14174
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 12 kDa; Observed MW: 12 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:50-1:100; FC-1:200-1:1000

Target information (P14174)

Gene Name	MIF
About Gene Symbol	



Western blot analysis of MIF in Jurkat lysates using MIF antibody.



Western blot analysis of MIF in C6, 3T3 lysates using MIF antibody.