

MMP2 Rabbit mAb


ABZ3055

Product Overview

Name	MMP2 Rabbit mAb
Accession(Primary)	P08253
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 74 kDa; Observed MW: 72,64 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100-1:200; FC-1:200-1:500

Target information (P08253)

Gene Name	MMP2
About Gene Symbol	



Western blot analysis of MMP2 in mouse lung,
mouse heart lysates using MMP2 antibody.