

MUC1 Rabbit mAb


ABZ3058

Product Overview

Name	MUC1 Rabbit mAb
Accession(Primary)	P15941
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 122 kDa; Observed MW: 25 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:200-1:500; FC-1:20; IP-1:20

Target information (P15941)

Gene Name	MUC1
About Gene Symbol	



Western blot analysis of MUC1 in HeLa, K562
lysates using Mucin 1 antibody.