

Muscarinic Acetylcholine Receptor M2 Rabbit mAb

ABZ3059


Product Overview

Name	Muscarinic Acetylcholine Receptor M2 Rabbit mAb
Accession(Primary)	P08172
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 52 kDa; Observed MW: 52 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200

Target information (P08172)

Gene Name CHRM2

About Gene Symbol



Western blot analysis of Muscarinic Acetylcholine Receptor 2/CM2 in K562, C6, HeLa lysates using Muscarinic Acetylcholine Receptor M2 antibody.

