


MOG Rabbit mAb

ABZ3060

Product Overview

Name	MOG Rabbit mAb
Accession(Primary)	Q16653
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 28 kDa; Observed MW: 28 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:200-1:500; IP-1:20-1:50

Target information (Q16653)

Gene Name MOG**About Gene Symbol**

Western blot analysis of Myelin oligodendrocyte glycoprotein in C6 lysates using MOG antibody.

