

Myeloperoxidase Rabbit mAb


ABZ3062

Product Overview


Name	Myeloperoxidase Rabbit mAb
Accession(Primary)	P05164
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 84 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:2000-1:20000; IHC-1:100-1:200; ICC/IF-1:50-1:100; FC-1:200-1:500

Target information (P05164)

Gene Name	MPO
About Gene Symbol	



Western blot analysis of Myeloperoxidase in rat spleen lysates using Myeloperoxidase antibody.



Western blot analysis of Myeloperoxidase in mouse spleen lysates using Myeloperoxidase antibody.