


# Myosin Light Chain 2 Rabbit mAb

ABZ3064

## Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | Myosin Light Chain 2 Rabbit mAb            |
| <b>Accession(Primary)</b>   | P10916                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                          |
| <b>Calculated MW</b>        | Calculated MW: 19 kDa; Observed MW: 19 kDa |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:200-1:500        |

## Target information (P10916)

**Gene Name** MYL2**About Gene Symbol**

Western blot analysis of Myosin Light Chain 2 in rat heart, mouse heart lysates using MYL2 antibody.

