

NAT1 Rabbit mAb


ABZ3065

Product Overview

Name	NAT1 Rabbit mAb
Accession(Primary)	P18440
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 34 kDa; Observed MW: 34 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100-1:200; FC-1:50-1:100

Target information (P18440)

Gene Name	NAT1
About Gene Symbol	



Western blot analysis of NAT1 in C6, 3T3, HeLa
lysates using NAT1 antibody.