

Nck Rabbit mAb


ABZ3067

Product Overview


Name	Nck Rabbit mAb
Accession(Primary)	P16333
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 43 kDa; Observed MW: 43 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:100-1:200; FC-1:200-1:1000

Target information (P16333)


Gene Name	NCK1
About Gene Symbol	



Immunocytochemistry analysis of NCK1 (green) in Hela using NCK1 antibody, and DAPI (blue).



Western blot analysis of Nck in K562, rat Brain, C6, 3T3, Hela lysates using Nck antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer tissue using Nck antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.