

NDUFB4 Rabbit mAb


ABZ3070

Product Overview


Name	NDUFB4 Rabbit mAb
Accession(Primary)	O95168
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 15 kDa; Observed MW: 15 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:20-1:50

Target information (O95168)


Gene Name	NDUFB4
About Gene Symbol	




Immunocytochemistry analysis of NDUFB4 (green) in HeLa using NDUFB4 antibody, and DAPI (blue).



Western blot analysis of NDUFB4 in mouse skeletal muscle lysates using NDUFB4 antibody.



Western blot analysis of NDUFB4 in K562, rat Brain lysates using NDUFB4 antibody



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using NDUFB4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.