


NDUFS6 Rabbit mAb

ABZ3073


Product Overview

Name	NDUFS6 Rabbit mAb
Accession(Primary)	O75380
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 14 kDa; Observed MW: 14 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200


Target information (O75380)

Gene Name NDUFS6**About Gene Symbol**

Western blot analysis of NDUFS6 in mouse heart lysates using NDUFS6 antibody.



Western blot analysis of NDUF56 in HeLa lysates using NDUF56 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using NDUF56 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.