


NEK4 Rabbit mAb

ABZ3075

Product Overview

Name	NEK4 Rabbit mAb
Accession(Primary)	P51957
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 95 kDa; Observed MW: 95 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:10-1:100

Target information (P51957)

Gene Name NEK4**About Gene Symbol**

Western blot analysis of NEK4 in Jurkat, HepG2, HeLa lysates using NEK4 antibody.

