

NELFE Rabbit mAb


ABZ3076

Product Overview


Name	NELFE Rabbit mAb
Accession(Primary)	P18615
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 43 kDa; Observed MW: 43 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:200-1:500; FC-1:200-1:500; IP-1:10-1:100; IHC-1:50-1:100

Target information (P18615)

Gene Name	NELFE
About Gene Symbol	



Western blot analysis of NELFe in C6 lysates using NELFE antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using NELFe antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.