

NeuN Rabbit mAb


ABZ3077

Product Overview


Name	NeuN Rabbit mAb
Accession(Primary)	A6NFN3
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 34 kDa; Observed MW: 46-55 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:500-1:2000; ICC/IF-1:100-1:200; FC-1:50-1:100

Target information (A6NFN3)


Gene Name	RBFOX3
About Gene Symbol	



Immunocytochemistry analysis of NeuN (green) in SH-SY5Y using NeuN antibody, and DAPI(blue).



Western blot analysis of NeuN in mouse brain lysates using NeuN antibody.



Western blot analysis of NeuN in rat Brain lysates
using NeuN antibody