

Neurofilament Heavy Polypeptide Rabbit mAb


ABZ3078

Product Overview

Name	Neurofilament Heavy Polypeptide Rabbit mAb
Accession(Primary)	P12036
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 112 kDa; Observed MW: 180-220 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:1000; ICC/IF-1:50-1:100

Target information (P12036)

Gene Name	NEFH
About Gene Symbol	



Western blot analysis of Neurofilament heavy polypeptide in mouse brain lysates using Neurofilament Heavy Polypeptide antibody.