

NFIB Rabbit mAb


ABZ3079

Product Overview

Name	NFIB Rabbit mAb
Accession(Primary)	O00712
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 47 kDa; Observed MW: 47,55 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1::50-1:100; FC-1:50-1:100; IP-1:20-1:50

Target information (O00712)

Gene Name	NFIB
About Gene Symbol	



Western blot analysis of NFIB in C6, HeLa lysates using NFIB antibody.