

NFYA Rabbit mAb


ABZ3080

Product Overview

Name	NFYA Rabbit mAb
Accession(Primary)	P23511
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 37 kDa; Observed MW: 37 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100-1:200; FC-1:10-1:100

Target information (P23511)

Gene Name	NFYA
About Gene Symbol	



Western blot analysis of NFYA in C6, 3T3, HeLa lysates using NFYA antibody.