

NR1D1 Rabbit mAb


ABZ3082

Product Overview


Name	NR1D1 Rabbit mAb
Accession(Primary)	P20393
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 67 kDa; Observed MW: 67 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:100-1:200; FC-1:20-1:50

Target information (P20393)

Gene Name	NR1D1
About Gene Symbol	



Western blot analysis of NR1D1 in K562, C6, 3T3, HeLa lysates using NR1D1 antibody.



Immunocytochemistry analysis of NR1D1(green) in HeLa using NR1D1 antibody, and DAPI(blue)