

NRG3 Rabbit mAb

ABZ3083


Product Overview

Name	NRG3 Rabbit mAb
Accession(Primary)	P56975
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 78 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:10-1:100

Target information (P56975)

Gene Name NRG3

About Gene Symbol



Western blot analysis of NRG3 in K562, rat Brain lysates using NRG3 antibody.

