

NSE Rabbit mAb


ABZ3084

Product Overview

Name	NSE Rabbit mAb
Accession(Primary)	P09104
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:50-1:100; FC-1:20; IP-1:20-1:50

Target information (P09104)

Gene Name	ENO2
About Gene Symbol	



Western blot analysis of NSE in Human Fetal Brain
lysates using NSE antibody.