

# Nucleophosmin Rabbit mAb


ABZ3085

## Product Overview


|                             |                                                                         |
|-----------------------------|-------------------------------------------------------------------------|
| <b>Name</b>                 | Nucleophosmin Rabbit mAb                                                |
| <b>Accession(Primary)</b>   | P06748                                                                  |
| <b>Host Species</b>         | Rabbit                                                                  |
| <b>Clonality</b>            | Monoclonal                                                              |
| <b>Isotype</b>              | IgG                                                                     |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                                       |
| <b>Calculated MW</b>        | Calculated MW: 33 kDa; Observed MW: 38 kDa                              |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:100-1:200;<br>ICC/IF-1:50-1:100; FC-1:20-1:50 |

## Target information (P06748)

|                          |      |
|--------------------------|------|
| <b>Gene Name</b>         | NPM1 |
| <b>About Gene Symbol</b> |      |



Western blot analysis of Nucleophosmin in HCT116, T47D, MCF-7 lysates using Nucleophosmin antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Nucleophosmin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.