

Nucleoside triphosphate phosphohydrolase Rabbit mAb


ABZ3086

Product Overview

Name	Nucleoside triphosphate phosphohydrolase Rabbit mAb
Accession(Primary)	Q9BSD7
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 21 kDa; Observed MW: 21 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:50-1:100; IP-1:20-1:50

Target information (Q9BSD7)

Gene Name	NTPCR
About Gene Symbol	



Western blot analysis of Nucleoside triphosphate
Phospho-hydrolase in K562, C6, 3T3, HeLa lysates
using NTPCR antibody.