


NUR77 Rabbit mAb

ABZ3090

Product Overview

Name	NUR77 Rabbit mAb
Accession(Primary)	P22736
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 64 kDa; Observed MW: 64 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P22736)

Gene Name NR4A1**About Gene Symbol**

Western blot analysis of NUR77 in rat Brain lysates using NUR77 antibody.

