

N WASP Rabbit mAb


ABZ3091

Product Overview


Name	N WASP Rabbit mAb
Accession(Primary)	O00401
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 55 kDa; Observed MW: 65 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:100-1:200; FC-1:100-1:200

Target information (O00401)


Gene Name	WASL
About Gene Symbol	



Western blot analysis of NWASP in 3T3 lysates using N WASP antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using NWASP antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of NWASP(green)
in Hela using NWASP antibody,and DAPI(blue)