

P Glycoprotein Rabbit mAb


ABZ3095

Product Overview

Name	P Glycoprotein Rabbit mAb
Accession(Primary)	P08183
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 141 kDa; Observed MW: 130-180 kDa
Recommended Dilution	WB-1:1000-1:50000; IHC-P-1:100-1:1200

Target information (P08183)

Gene Name	ABCB1
About Gene Symbol	ABCB1 (ATP Binding Cassette Subfamily B Member 1) is a cell surface marker involved in adaptive immunity. It is an important research target in immunology, flow cytometry, and cancer immunotherapy and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of P Glycoprotein in HeLa lysates using P Glycoprotein antibody.