

p40 Rabbit mAb


ABZ3096

Product Overview

Name	p40 Rabbit mAb
Accession(Primary)	Q7Z6M1
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 41 kDa; Observed MW: 41 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q7Z6M1)

Gene Name	RABEPK
About Gene Symbol	RABEPK (Rab9 Effector Protein With Kelch Motifs) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of p40 in K562, HeLa lysates using p40 antibody.