

P4HB Rabbit mAb


ABZ3098

Product Overview


Name	P4HB Rabbit mAb
Accession(Primary)	P07237
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 57 kDa; Observed MW: 57 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100

Target information (P07237)

Gene Name	P4HB
About Gene Symbol	P4HB (Prolyl 4-Hydroxylase Subunit Beta) is a collagen extracellular matrix protein involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular



Western blot analysis of P4HB in K562, C6 lysates using P4HB antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using P4HB antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.