

p53 Rabbit mAb


ABZ3099

Product Overview

Name	p53 Rabbit mAb
Accession(Primary)	P02340
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse
Calculated MW	Calculated MW: 43 kDa; Observed MW: 53 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:50-1:100; FC-1:50-1:100; IP-1:20-1:50

Target information (P02340)

Gene Name	Tp53
About Gene Symbol	



Western blot analysis of p53 in 3T3 lysates using p53 antibody.