

p60 Katanin Rabbit mAb


ABZ3100

Product Overview

Name	p60 Katanin Rabbit mAb
Accession(Primary)	O75449
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 56 kDa; Observed MW: 56 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100

Target information (O75449)

Gene Name	KATNA1
About Gene Symbol	KATNA1 (Katanin Catalytic Subunit A1) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of p60 katanin in K562, HeLa lysates using p60 Katanin antibody.