

p73 Rabbit mAb


ABZ3101

Product Overview

Name	p73 Rabbit mAb
Accession(Primary)	O15350
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 70 kDa; Observed MW: 80,70 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20-1:50

Target information (O15350)

Gene Name	TP73
About Gene Symbol	TP73 (Tumor Protein P73) is a transcription factor involved in gene transcription and regulation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of p73 in 293 lysates using p73 antibody.