

Nibrin Rabbit mAb


ABZ3102

Product Overview


Name	Nibrin Rabbit mAb
Accession(Primary)	O60934
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 85 kDa; Observed MW: 95 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:50-1:100

Target information (O60934)

Gene Name	NBN
About Gene Symbol	



Western blot analysis of p95/NBS1 in HeLa lysates using Nibrin antibody.



Immunocytochemistry analysis of p95/NBS1 (green) in MCF-7 using p95/NBS1 antibody, and DAPI (blue)