

PRKRA Rabbit mAb


ABZ3103

Product Overview

Name	PRKRA Rabbit mAb
Accession(Primary)	O75569
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 34 kDa; Observed MW: 34 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:100-1:200

Target information (O75569)

Gene Name	PRKRA
About Gene Symbol	PRKRA (Protein Activator Of Interferon Induced Protein Kinase EIF2AK2) is a cytokine and immune signaling protein involved in JAK-STAT/cytokine signaling and innate immune signaling. It is an important research target in immunology, inflammation, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PACT (Phospho- KR activating protein)/PRKRA in C6, 3T3, HeLa lysates using PRKRA antibody.