

PAI1 Rabbit mAb


ABZ3104

Product Overview

Name	PAI1 Rabbit mAb
Accession(Primary)	P05121
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 45 kDa; Observed MW: 45 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:50-1:100; FC-1:50-1:100; IP-1:20-1:50

Target information (P05121)

Gene Name	SERPINE1
About Gene Symbol	SERPINE1 (Serpin Family E Member 1) is a protease enzyme involved in angiogenesis. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PAI1 in 3T3, HeLa lysates using PAI1 antibody.