

P cadherin Rabbit mAb


ABZ3106

Product Overview


Name	P cadherin Rabbit mAb
Accession(Primary)	P22223
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 97 kDa; Observed MW: 130-150 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:100-1:200; FC-1:20-1:50; IP-1:20

Target information (P22223)


Gene Name	CDH3
About Gene Symbol	




Western blot analysis of pan Cadherin in mouse heart, mouse brain lysates using Cadherin antibody.



Western blot analysis of pan Cadherin in HeLa, CHO-K1, C6, rat Brain lysates using pan Cadherin antibody.



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using pan Cadherin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of pan Cadherin (green) in MCF-7 using pan Cadherin antibody, and DAPI (blue).

