

PARP1 Rabbit mAb


ABZ3107

Product Overview


Name	PARP1 Rabbit mAb
Accession(Primary)	P09874
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 113 kDa; Observed MW: 116,89 kDa
Recommended Dilution	WB-1:1000-1:10000; IHC-1:100-1:1000; ICC/IF-1:200-1:500

Target information (P09874)


Gene Name	PARP1
About Gene Symbol	



Western blot analysis of PARP1 in K562, 3T3, HeLa lysates using PARP1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using PARP1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of PARP1 (green) in HeLa using PARP1 antibody, and DAPI (blue)

