


PCPTP1 Rabbit mAb

ABZ3108

Product Overview

Name	PCPTP1 Rabbit mAb
Accession(Primary)	Q15256
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 74 kDa; Observed MW: 74 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q15256)

Gene Name PTPRR**About Gene Symbol**

Western blot analysis of PCPTP1/PTPRQ in rat
Brain lysates using PCPTP1 antibody.

