

PDGFR alpha Rabbit mAb


ABZ3111

Product Overview

Name	PDGFR alpha Rabbit mAb
Accession(Primary)	P16234
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 123 kDa; Observed MW: 190 kDa
Recommended Dilution	WB-1:1000-1:10000; IHC-1:50-1:200; ICC/IF-1:200-1:500; FC-1:20-1:50; IP-1:20

Target information (P16234)

Gene Name	PDGFRA
About Gene Symbol	



Western blot analysis of PDGF Receptor alpha in 3T3 lysates using PDGFR alpha antibody.