

PDLIM1/CLP36 Rabbit mAb


ABZ3112

Product Overview


Name	PDLIM1/CLP36 Rabbit mAb
Accession(Primary)	O00151
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 36 kDa; Observed MW: 36 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:50-1:100; FC-1:100-1:500

Target information (O00151)

Gene Name	PDLIM1
About Gene Symbol	PDLIM1 (PDZ And LIM Domain 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PDLIM1 in K562, C6, 3T3, HeLa lysates using PDLIM1 antibody.



Immunocytochemistry analysis of PDLIM1(green) in HeLa using PDLIM1 antibody, and DAPI(blue)