

PEF1 Rabbit mAb


ABZ3113

Product Overview

Name	PEF1 Rabbit mAb
Accession(Primary)	Q9UBV8
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 30 kDa; Observed MW: 30 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:10-1:100

Target information (Q9UBV8)

Gene Name	PEF1
About Gene Symbol	PEF1 (Penta-EF-Hand Domain Containing 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PEF1 in K562, rat Brain, C6, 3T3, HeLa lysates using PEF1 antibody.