


PEX11B Rabbit mAb

ABZ3115


Product Overview

Name	PEX11B Rabbit mAb
Accession(Primary)	O96011
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 28 kDa; Observed MW: 28 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (O96011)

Gene Name PEX11B**About Gene Symbol**

Western blot analysis of PEX11B in T47D, HL-60 lysates using PEX11B antibody.



Western blot analysis of PEX11B in 3T3 lysates
using PEX11B antibody.