

UCH-L1/PGP9.5 Rabbit mAb


ABZ3116

Product Overview


Name	UCH-L1/PGP9.5 Rabbit mAb
Accession(Primary)	P09936
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 25 kDa; Observed MW: 25 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:200-1:500; FC-1:50-1:100; IP-1:10-1:100

Target information (P09936)

Gene Name	UCHL1
About Gene Symbol	UCHL1 (Ubiquitin C-Terminal Hydrolase L1) is a protein kinase and signaling regulator involved in ubiquitin-proteasome pathway. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PGP9.5 in 293, SH-SY5Y, mouse brain lysates using PGP9.5 antibody.



Western blot analysis of PGP9.5 in rat Brain lysates using PGP9.5 antibody.