

PIN4 Rabbit mAb


ABZ3117

Product Overview

Name	PIN4 Rabbit mAb
Accession(Primary)	Q9Y237
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 14 kDa; Observed MW: 14 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; IF-1:50-1:100

Target information (Q9Y237)

Gene Name	PIN4
About Gene Symbol	



Western blot analysis of PIN4 in CHO-K1, C6, Jurkat lysates using PIN4 antibody.