

PIP5K1C Rabbit mAb

ABZ3118


Product Overview

Name	PIP5K1C Rabbit mAb
Accession(Primary)	O70161
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 72 kDa; Observed MW: 90,87 kDa
Recommended Dilution	WB-1:1000-1:5000


Target information (O70161)

Gene Name Pip5K1C

About Gene Symbol



Western blot analysis of PIP5K1 gamma in mouse brain lysates using PIP5K1C antibody.



Western blot analysis of PIP5KI gamma in MCF-7
lysates using PIP5KI gamma antibody.