

PKC beta 2 Rabbit mAb


ABZ3119

Product Overview



Name	PKC beta 2 Rabbit mAb
Accession(Primary)	P05771
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 77 kDa; Observed MW: 77 kDa
Recommended Dilution	WB-1:500-1:1000; IHC:1:20; ICC/IF-1:20-1:50; FC-1:20; IP-1:20

Target information (P05771)


Gene Name	PRKCB
About Gene Symbol	



Immunocytochemistry analysis of PKC beta 2 (green) in K562 using PKC beta 2 antibody, and DAPI(blue).



Western blot analysis of PKC beta 2 in rat brain lysates using PKC beta 2 antibody.



Western blot analysis of PKC beta 2 in K562
lysates using PKC beta 2 antibody